

For Research Use Only

NOP10 Recombinant monoclonal antibody

Catalog Number:87529-1-RR



Basic Information

Catalog Number: 87529-1-RR	GenBank Accession Number: BC063023	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 55505	CloneNo.: 252910G7
Isotype: IgG	UNIPROT ID: Q9NPE3	Recommended Dilutions: WB: 1:5000-1:50000 IHC: 1:50-1:500
Immunogen Catalog Number: AG18718	Full Name: NOP10 ribonucleoprotein homolog (yeast)	
	Calculated MW: 64 aa, 77 kDa	
	Observed MW: 10-12 kDa	

Applications

Tested Applications:

WB, IHC, ELISA

Species Specificity:

human, mouse, rat

Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0

Positive Controls:

WB : HeLa cells, Jurkat cells, Caco-2 cells, COLO 320 cells, mouse thymus tissue, rat thymus tissue

IHC : Human brain cancer (Astrocytoma) ,

Background Information

Nucleolar protein 10 (NOP10) is a component of the H/ACA small nucleolar ribonucleoprotein (snoRNP) complex and is essential for maintaining telomerase activity and RNA pseudouridylation. NOP10 is abnormally expressed in a variety of cancers such as breast cancer, pancreatic cancer, lung cancer, etc., and high expression is usually associated with poor prognosis (PMID: 33161513; 39538226).

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 87529-1-RR (NOP10 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.



Immunohistochemical analysis of paraffin-embedded Human brain cancer (Astrocytoma) slide using 87529-1-RR (NOP10 antibody) at dilution of 1:200 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunohistochemical analysis of paraffin-embedded Human brain cancer (Astrocytoma) slide using 87529-1-RR (NOP10 antibody) at dilution of 1:200 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).