

For Research Use Only

TEKT1 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 87530-2-PBS



Basic Information

Catalog Number:

87530-2-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG5173

GenBank Accession Number:

BC014599

GeneID (NCBI):

83659

UNIPROT ID:

Q969V4

Full Name:

tektin 1

Calculated MW:

48 kDa

Purification Method:

Protein A purification

CloneNo.:

252780D10

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean TEK1 concentration was determined to be 1691.4 ng/mL in mouse testis tissue extract based on a 4.2 mg/mL extract load, mouse pancreas tissue extract were served as a negative control.



Cytometric bead array standard curve of MP03082-1, TEK1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87530-3-PBS. Detection antibody: 87530-2-PBS. Standard: Ag5173. Range: 0.781-100 ng/mL