

For Research Use Only

TEKT1 Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number:87530-3-PBS



Basic Information

Catalog Number: 87530-3-PBS	GenBank Accession Number: BC014599	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 83659	CloneNo.: 252780E10
Isotype: IgG	UNIPROT ID: Q969V4	
Immunogen Catalog Number: AG5173	Full Name: tektin 1	
	Calculated MW: 48 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean TEK1 concentration was determined to be 1691.4 ng/mL in mouse testis tissue extract based on a 4.2 mg/mL extract load, mouse pancreas tissue extract were served as a negative control.



Cytometric bead array standard curve of MP03082-1, TEK1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87530-3-PBS. Detection antibody: 87530-2-PBS. Standard: Ag5173. Range: 0.781-100 ng/mL