

For Research Use Only

TBP Recombinant monoclonal antibody, PBS Only

Catalog Number: 87550-1-PBS



Basic Information

Catalog Number: 87550-1-PBS	GenBank Accession Number: BC110341	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 6908	CloneNo.: 252942E3
Isotype: IgG	UNIPROT ID: P20226	
Immunogen Catalog Number: AG12383	Full Name: TATA box binding protein	
	Calculated MW: 338 aa, 38 kDa	
	Observed MW: 37-45 kDa	

Applications

Tested Applications:
WB, IF/ICC, Indirect ELISA

Species Specificity:
human

Background Information

TBP, TATA box-binding protein, is a DNA-binding subunit of TFIID, which is a multisubunit complex essential for the expression of most protein-encoding genes. The well-conserved C-terminal domains of 180 aa is both necessary and sufficient for interaction with DNA and for assembly of the basal transcription apparatus. The N terminus of TBP modulates the DNA-binding activity of the C terminus of the protein. Once binds to the TATA box mediated by TBP, TFIID initiate the transcription step of the pre-initiation complex (PIC), playing a role in the activation of eukaryotic genes transcribed by RNA polymerase II.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 87550-1-RR (TBP antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 87550-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (4% PFA) fixed HepG2 cells using TBP antibody (87550-1-RR, Clone: 252942E3) at dilution of 1:1000 and CoraLite® 488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2), CL594-Phalloidin (red). This data was developed using the same antibody clone with 87550-1-PBS in a different storage buffer formulation.