

For Research Use Only

Rat FABP1 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 87555-1-PBS



Basic Information

Catalog Number: 87555-1-PBS	GenBank Accession Number: NM_012556.2	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 24360	CloneNo.: 252724F10
Isotype: IgG	UNIPROT ID: P02692	
Immunogen Catalog Number: EG7277	Full Name: fatty acid binding protein 1, liver	
	Calculated MW: 14 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
rat

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean FABP1 concentration was determined to be 44.4 ng/mL in rat serum.



Cytometric bead array standard curve of MP03047-1, RAT FABP1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87555-2-PBS. Detection antibody: 87555-1-PBS. Standard: Eg7277. Range: 0.156-20 ng/mL.