

For Research Use Only

PTPN23 Recombinant monoclonal antibody, PBS Only

Catalog Number: 87561-1-PBS



Basic Information

Catalog Number: 87561-1-PBS	GenBank Accession Number: BC004881	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 25930	CloneNo.: 252957F3
Isotype: IgG	UNIPROT ID: Q9H3S7	
Immunogen Catalog Number: AG41272	Full Name: protein tyrosine phosphatase, non-receptor type 23	
	Calculated MW: 179 kDa	
	Observed MW: 179-200 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 87561-1-RR (PTPN23 antibody) at dilution of 1:20000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 87561-1-PBS in a different storage buffer formulation.