

For Research Use Only

GFRA3 Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 87564-2-PBS



Basic Information

Catalog Number: 87564-2-PBS	GenBank Accession Number: NM_001496	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 2676	CloneNo.: 252618A3
Isotype: IgG	UNIPROT ID: O60609-1	
Immunogen Catalog Number: EG5897	Full Name: GDNF family receptor alpha 3	
	Calculated MW: 45 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean GFRA3 concentration was determined to be 41.0 ng/mL in SH-SY5Y cell extract based on a 1.2 mg/mL extract load, HEK-293 cell extract were served as a negative control.



Cytometric bead array standard curve of MP03049-1, GFRA3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87564-2-PBS. Detection antibody: 87564-1-PBS. Standard: Eg5897. Range: 0.312-40 ng/mL.