

For Research Use Only

CNTN2 Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 87574-2-PBS



Basic Information

Catalog Number: 87574-2-PBS	GenBank Accession Number: NM_001346083.2	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 6900	CloneNo.: 252726B2
Isotype: IgG	UNIPROT ID: Q02246	
Immunogen Catalog Number: EG5489	Full Name: contactin 2 (axonal)	
	Calculated MW: 113 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean CNTN2 concentration was determined to be 168.6 ng/mL in mouse brain tissue extract based on a 3.2 mg/mL extract load, mouse liver tissue extract were served as a negative control.



Cytometric bead array standard curve of MP03102-1, CNTN2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87574-2-PBS. Detection antibody: 87574-1-PBS. Standard: Eg5489. Range: 0.625-80 ng/mL.