

For Research Use Only

CD296/ART1 Recombinant monoclonal antibody

Catalog Number: 87577-1-RR



Basic Information

Catalog Number: 87577-1-RR	GenBank Accession Number: NM_009710.4	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 11870	CloneNo.: 252922A3
Isotype: IgG	UNIPROT ID: Q60935	Recommended Dilutions: WB: 1:5000-1:50000
Immunogen Catalog Number: EG6314	Full Name: ADP-ribosyltransferase 1	
	Calculated MW: 36 kDa	
	Observed MW: 40 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls:
Species Specificity: mouse, rat	WB : mouse skeletal muscle tissue, rat skeletal muscle tissue, mouse heart tissue

Background Information

CD296 (ART1) is a GPI-anchored ecto-ADP-ribosyltransferase that catalyzes extracellular mono-ADP-ribosylation of cell-surface proteins. By modifying integrins and other membrane-associated targets, ART1 regulates cell adhesion, migration, and inflammatory signaling. Its overexpression in multiple cancers and involvement in airway and tissue remodeling highlight ART1 as a functional mediator at the interface of cell-surface modification and disease progression.

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 87577-1-RR (CD296/ART1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.



Biolayer interferometry (BLI) kinetic assays of 87577-1-RR against Mouse CD296/ART1 were performed. The affinity constant is 4.56 nM.