

For Research Use Only

Mouse CD296/ART1 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 87577-2-PBS



Basic Information

Catalog Number: 87577-2-PBS	GenBank Accession Number: NM_009710.4	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 11870	CloneNo.: 252922F2
Isotype: IgG	UNIPROT ID: Q60935	
Immunogen Catalog Number: EG6314	Full Name: ADP-ribosyltransferase 1	
	Calculated MW: 36 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
mouse

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean ART1 concentration was determined to be 37.6 ng/mL in mouse heart tissue extract based on a 5.4 mg/mL extract load, mouse brain tissue extract were served as a negative control.



Cytometric bead array standard curve of MP03128-1, MOUSE ART1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87577-3-PBS. Detection antibody: 87577-2-PBS. Standard: Eg6314. Range: 0.078-10 ng/mL.