

For Research Use Only

TMEM150 Recombinant monoclonal antibody, PBS Only

Catalog Number: 87578-1-PBS



Basic Information

| | | |
|---|--|---|
| Catalog Number: 87578-1-PBS | GenBank Accession Number: BC050466 | Purification Method: Protein A purification |
| Source: Rabbit | GeneID (NCBI): 129303 | CloneNo.: 252959H6 |
| Isotype: IgG | UNIPROT ID: Q86TG1 | |
| Immunogen Catalog Number: AG36287 | Full Name: transmembrane protein 150 | |
| | Calculated MW: 271 aa, 29 kDa | |
| | Observed MW: 29 kDa | |

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human, mouse, rat

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 87578-1-RR (TMEM150 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 87578-1-PBS in a different storage buffer formulation.