

For Research Use Only

# TMEM150 Recombinant monoclonal antibody

Catalog Number:87578-1-RR



## Basic Information

Catalog Number:

87578-1-RR

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG36287

GenBank Accession Number:

BC050466

GeneID (NCBI):

129303

UNIPROT ID:

Q86TG1

Full Name:

transmembrane protein 150

Calculated MW:

271 aa, 29 kDa

Observed MW:

29 kDa

Purification Method:

Protein A purification

CloneNo.:

252959H6

Recommended Dilutions:

WB: 1:500-1:2000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : Jurkat cells, Neuro-2a cells, PC-12 cells

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 87578-1-RR (TMEM150 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.