

For Research Use Only

CBFA2T2 Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 87598-4-PBS



Basic Information

Catalog Number: 87598-4-PBS	GenBank Accession Number: BC015066	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 9139	CloneNo.: 252968G11
Isotype: IgG	UNIPROT ID: O43439	
Immunogen Catalog Number: AG33484	Full Name: core-binding factor, runt domain, alpha subunit 2; translocated to, 2	
	Calculated MW: 604 aa, 67 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean CBFA2T2 concentration was determined to be 46.6 ng/mL in SH-SY5Y cell extract based on a 1.2 mg/mL extract load, 32.3 ng/mL in HEK-293T cell extract based on a 1.2 mg/mL extract load, 3.8 ng/mL in HeLa cell extract based on a 1.2 mg/mL extract load.



Cytometric bead array standard curve of MP03122-1, CBFA2T2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87598-4-PBS. Detection antibody: 87598-3-PBS. Standard: Ag33484. Range: 0.781-100 ng/mL