

For Research Use Only

# PEG Recombinant monoclonal antibody

Catalog Number:87604-1-RR



## Basic Information

Catalog Number:

87604-1-RR

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

GeneID (NCBI):

Full Name:

Purification Method:

Protein A purification

CloneNo.:

240925A1

Recommended Dilutions:

ELISA: 1:10-1:100

## Applications

Tested Applications:

ELISA

Species Specificity:

-

Positive Controls:

ELISA : PEG,

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



96-well plates were coated with Peptid PEG12, mPEG20K-Biotin, Biotin-PEG2K-Biotin. The plate was blocked with 1% BSA blocking buffer, then washed, and incubated with 100  $\mu$ L per well of PEG antibody 87604-1-RR and serially diluted with a starting concentration of 2  $\mu$ g/mL (1 hr, 37°C). The plate was washed and incubated with 100  $\mu$ L per well of an HRP-conjugated goat anti-Rabbit IgG secondary antibody at a dilution of 1:10,000 (1 hr, 37°C). Detection was performed using an Ultra TMB

