

For Research Use Only

Rat CXCL16 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number:87605-1-PBS



Basic Information

Catalog Number: 87605-1-PBS	GenBank Accession Number: NM_001017478.1	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 497942	CloneNo.: 252531D1
Isotype: IgG	UNIPROT ID: Q6AXU5	
Immunogen Catalog Number: EG6251	Full Name: chemokine (C-X-C motif) ligand 16	
	Calculated MW: 27 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
rat

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean CXCL16 concentration was determined to be 5.9 ng/mL in rat serum.



Cytometric bead array standard curve of MP03078-1, RAT CXCL16 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87605-2-PBS. Detection antibody: 87605-1-PBS. Standard: Eg6251. Range: 0.078-10 ng/mL.