

For Research Use Only

Human IgA1 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number:87606-1-PBS



Basic Information

Catalog Number:

87606-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG9190

GenBank Accession Number:

GeneID (NCBI):

3493

Full Name:

immunoglobulin heavy constant
alpha 1

Purification Method:

Protein A purification

CloneNo.:

252592G2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean IgA1 concentration was determined to be 2339.4 ug/mL in human serum.



Cytometric bead array standard curve of MP03080-1, IgA1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87606-2-PBS. Detection antibody: 87606-1-PBS. Standard: Eg9190. Range: 0.156-20 ng/mL.