

For Research Use Only

Mouse Abcg2 Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number:87614-3-PBS



Basic Information

Catalog Number:

87614-3-PBS

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC053730

GeneID (NCBI):

26357

UNIPROT ID:

Q7TMS5

Full Name:

ATP-binding cassette, sub-family G (WHITE), member 2

Purification Method:

Protein A purification

CloneNo.:

252770C9

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean BCRP/ABCG2 concentration was determined to be 8.7 ng/mL in mouse liver tissue extract based on a 8.4 mg/mL extract load, mouse eye tissue extract were served as a negative control.



Cytometric bead array standard curve of MP03130-1, MOUSE BCRP/ABCG2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87614-3-PBS. Detection antibody: 87614-2-PBS. Standard: SY02955. Range: 0.156-20 ng/mL