

For Research Use Only

ZBTB2 Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number:87621-2-PBS



Basic Information

Catalog Number: 87621-2-PBS	GenBank Accession Number: BC020172	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 57621	CloneNo.: 252923C6
Isotype: IgG	UNIPROT ID: Q8N680	
Immunogen Catalog Number: AG24041	Full Name: zinc finger and BTB domain containing 2	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean ZBTB2 concentration was determined to be 10.3 ng/mL in Raji cell extract based on a 1.2 mg/mL extract load, 3.8 ng/mL in MCF-7 cell extract based on a 1.2 mg/mL extract load, 1.9 ng/mL in COLO 320 cell extract based on a 1.2 mg/mL extract load.



Cytometric bead array standard curve of MP03125-1, ZBTB2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87621-2-PBS. Detection antibody: 87621-1-PBS. Standard: Ag24041. Range: 0.781-100 ng/mL