

For Research Use Only

ACO2 Recombinant monoclonal antibody, PBS Only

Catalog Number:87630-2-PBS



Basic Information

Catalog Number: 87630-2-PBS	GenBank Accession Number: BC014092	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 50	CloneNo.: 252976A4
Isotype: IgG	UNIPROT ID: Q99798	
Immunogen Catalog Number: AG1590	Full Name: aconitase 2, mitochondrial	
	Calculated MW: 85 kDa	
	Observed MW: 85 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human

Background Information

ACO2(aconitate hydratase, mitochondrial) is also named as citrate hydro-lyase and belongs to the aconitase/IPM isomerase family. It plays a key function in cellular energy production, and loss of its activity has a major impact on cellular and organismal survival.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 87630-2-RR (ACO2 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 87630-2-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 87630-2-RR against Human ACO2 were performed. The affinity constant is below 1 pM.