

For Research Use Only

TEKT2 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number:87634-2-PBS



Basic Information

Catalog Number: 87634-2-PBS	GenBank Accession Number: BC035620	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 27285	CloneNo.: 252977G12
Isotype: IgG	UNIPROT ID: Q9UIF3	
Immunogen Catalog Number: AG4457	Full Name: tektin 2 (testicular)	
	Calculated MW: 430 aa, 50 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean TEKT2 concentration was determined to be 49.2 ng/mL in mouse testis tissue extract based on a 4.2 mg/mL extract load, 2.1 ng/mL in mouse lung tissue extract based on a 5.0 mg/mL extract load.



Cytometric bead array standard curve of MP03140-1, TEKT2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87634-3-PBS. Detection antibody: 87634-2-PBS. Standard: Ag4457. Range: 0.625-80 ng/mL.