

For Research Use Only

Mouse TLT-2/Treml2 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 87644-1-PBS



Basic Information

Catalog Number: 87644-1-PBS	GenBank Accession Number: NM_001033405	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 328833	CloneNo.: 252927C7
Isotype: IgG	UNIPROT ID: Q2LA85	
Immunogen Catalog Number: EG7510	Full Name: triggering receptor expressed on myeloid cells-like 2	
	Calculated MW: 37 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
mouse

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean TLT-2/Trem12 concentration was determined to be 90.8 ng/mL in mouse spleen tissue extract based on a 6.4 mg/mL extract load, mouse heart tissue extract were served as a negative control.



Cytometric bead array standard curve of MP03105-1, MOUSE TLT-2/Trem12 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87644-2-PBS. Detection antibody: 87644-1-PBS. Standard: Eg7510. Range: 0.469-60 ng/mL.