

For Research Use Only

ZNF207 Recombinant monoclonal antibody, PBS Only

Catalog Number:87649-1-PBS



Basic Information

| | | |
|---|--|---|
| Catalog Number: 87649-1-PBS | GenBank Accession Number: BC002372 | Purification Method: Protein A purification |
| Source: Rabbit | GeneID (NCBI): 7756 | CloneNo.: 253049C10 |
| Isotype: IgG | UNIPROT ID: O43670 | |
| Immunogen Catalog Number: AG40360 | Full Name: zinc finger protein 207 | |
| | Calculated MW: 51 kDa | |
| | Observed MW: 50-57 kDa | |

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human, mouse, rat

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

