

For Research Use Only

ZNF207 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 87649-2-PBS



Basic Information

Catalog Number: 87649-2-PBS	GenBank Accession Number: BC002372	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 7756	CloneNo.: 253049B3
Isotype: IgG	UNIPROT ID: O43670	
Immunogen Catalog Number: AG40360	Full Name: zinc finger protein 207	
	Calculated MW: 51 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean ZNF207 concentration was determined to be 6.2 ng/mL in Ramos cell extract based on a 1.2 mg/mL extract load, 1.5 ng/mL in PC-3 cell extract based on a 1.2 mg/mL extract load.



Cytometric bead array standard curve of MP03134-1, ZNF207 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87649-3-PBS. Detection antibody: 87649-2-PBS. Standard: Ag40360. Range: 0.078-10 ng/mL.