

For Research Use Only

# ZNF207 Recombinant monoclonal antibody, PBS Only (Capture/Detector)

Catalog Number: 87649-3-PBS



## Basic Information

<b>Catalog Number:</b> 87649-3-PBS	<b>GenBank Accession Number:</b> BC002372	<b>Purification Method:</b> Protein A purification
<b>Source:</b> Rabbit	<b>GeneID (NCBI):</b> 7756	<b>CloneNo.:</b> 253049E6
<b>Isotype:</b> IgG	<b>UNIPROT ID:</b> O43670	
<b>Immunogen Catalog Number:</b> AG40360	<b>Full Name:</b> zinc finger protein 207	
	<b>Calculated MW:</b> 51 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Sandwich ELISA, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**


## Selected Validation Data




The mean ZNF207 concentration was determined to be 8.0 ng/mL in Ramos cell extract based on a 1.5 mg/mL extract load, 4.6 ng/mL in Raji cell extract based on a 1.2 mg/mL extract load and 8.2 ng/mL in PC-3 cell extract based on a 1.5 mg/mL extract load.



Sandwich ELISA standard curve of MP03134-2, Human ZNF207 Recombinant Matched Antibody Pair - PBS only. 87649-4-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag40360. 87649-3-PBS was HRP conjugated as the detection antibody. Range: 78.1-5000 pg/mL



The mean ZNF207 concentration was determined to be 6.2 ng/mL in Ramos cell extract based on a 1.2 mg/mL extract load, 1.5 ng/mL in PC-3 cell extract based on a 1.2 mg/mL extract load.



Cytometric bead array standard curve of MP03134-1, ZNF207 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87649-3-PBS. Detection antibody: 87649-2-PBS. Standard: Ag40360. Range: 0.078-10 ng/mL