

For Research Use Only

# Periostin Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number:87652-1-PBS



## Basic Information

Catalog Number:

87652-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG4564

GenBank Accession Number:

NM\_006475.3

GeneID (NCBI):

10631

UNIPROT ID:

Q15063-1

Full Name:

periostin, osteoblast specific factor

Calculated MW:

93 kDa

Purification Method:

Protein A purification

CloneNo.:

252916G3

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



The mean Periostin concentration was determined to be 1226.7 ng/mL in human serum.



Cytometric bead array standard curve of MP03123-1, Periostin Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87652-2-PBS. Detection antibody: 87652-1-PBS. Standard: Eg4564. Range: 1.562-200 ng/mL