

For Research Use Only

# RAB3D Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 87655-3-PBS



## Basic Information

Catalog Number:

87655-3-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG5570

GenBank Accession Number:

NM\_004283.3

GeneID (NCBI):

9545

UNIPROT ID:

O95716

Full Name:

RAB3D, member RAS oncogene family

Calculated MW:

24 kDa

Purification Method:

Protein A purification

CloneNo.:

253078D8

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



The mean RAB3D concentration was determined to be 15.4 ng/mL in SK-BR-3 cell extract based on a 1.5 mg/mL extract load, HT-1080 cell extract were served as a negative control.



Cytometric bead array standard curve of MP03138-1, RAB3D Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87655-2-PBS. Detection antibody: 87655-3-PBS. Standard: Eg5570. Range: 0.156-20 ng/mL.