

For Research Use Only

# Mouse Ccl28 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number:87665-1-PBS



## Basic Information

Catalog Number:

87665-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG5201

GenBank Accession Number:

NM\_020279.3

GeneID (NCBI):

56838

UNIPROT ID:

Q9JIL2

Full Name:

chemokine (C-C motif) ligand 28

Calculated MW:

15 kDa

Purification Method:

Protein A purification

CloneNo.:

252525D2

## Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

mouse

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



The mean CCL28 concentration was determined to be 2,422.24 ng/mL in mouse testis tissue extract based on a 4.2 mg/mL extract load and 1,834.39 ng/mL in mouse colon tissue extract based on a 4.5 mg/mL extract load.



Sandwich ELISA standard curve of MP03113-1, Mouse CCL28 Recombinant Matched Antibody Pair - PBS only. 87665-2-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg5201. 87665-1-PBS was HRP conjugated as the detection antibody. Range: 0.78-50 ng/mL