

For Research Use Only

Mouse Ccl28 Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number:87665-2-PBS



Basic Information

Catalog Number: 87665-2-PBS	GenBank Accession Number: NM_020279.3	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 56838	CloneNo.: 252525A7
Isotype: IgG	UNIPROT ID: Q9JIL2	
Immunogen Catalog Number: EG5201	Full Name: chemokine (C-C motif) ligand 28	
	Calculated MW: 15 kDa	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA

Species Specificity:
mouse

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean CCL28 concentration was determined to be 2,422.24 ng/mL in mouse testis tissue extract based on a 4.2 mg/mL extract load and 1,834.39 ng/mL in mouse colon tissue extract based on a 4.5 mg/mL extract load.



Sandwich ELISA standard curve of MP03113-1, Mouse CCL28 Recombinant Matched Antibody Pair - PBS only. 87665-2-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg5201. 87665-1-PBS was HRP conjugated as the detection antibody. Range: 0.78-50 ng/mL