

For Research Use Only

PAX9 Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number:87678-2-PBS



Basic Information

Catalog Number:

87678-2-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG25381

GenBank Accession Number:

BC001159

GeneID (NCBI):

5083

UNIPROT ID:

P55771

Full Name:

paired box 9

Calculated MW:

341 aa, 36 kDa

Purification Method:

Protein A purification

CloneNo.:

253236E9

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean PAX9 concentration was determined to be 9.6 ng/mL in A2780 cell extract based on a 1.2 mg/mL extract load, SH-SY5Y cell extract were served as a negative control.



Cytometric bead array standard curve of MP03150-1, PAX9 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87678-2-PBS. Detection antibody: 87678-1-PBS. Standard: Ag25381. Range: 0.156-20 ng/mL.