

For Research Use Only

Collagen Type XXIII Recombinant monoclonal antibody, PBS Only

Catalog Number:87687-3-PBS



Basic Information

Catalog Number: 87687-3-PBS	GenBank Accession Number: NM_173465.3	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 91522	CloneNo.: 253143A9
Isotype: IgG	UNIPROT ID: Q86Y22	
Immunogen Catalog Number: EG4452	Full Name: collagen, type XXIII, alpha 1	
	Calculated MW: 52 kDa	
	Observed MW: 68-75 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human, mouse, rat

Background Information

Collagen XXIII is a type II transmembrane protein and the fourth in the subfamily of non-fibrillar transmembranous collagens. This kind of collagens have a single pass hydrophobic transmembrane domain. The molecule of collagen XXIII can be found either in membrane-bound form or in shed form. Type XXIII collagen is expressed in both adult tissues and developing organs. It can be found in the epidermis and other epithelia such as those in tongue, gut and lung, but also in the brain, the kidney and the cornea. It has been shown that in prostate collagen XXIII expression is associated with tumor progression.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 87687-3-RR (Collagen XXIII alpha 1/ COL23A1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 87687-3-PBS in a different storage buffer formulation.