

For Research Use Only

# MARCS Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number:87693-2-PBS



## Basic Information

<b>Catalog Number:</b> 87693-2-PBS	<b>GenBank Accession Number:</b> NM_002356.7	<b>Purification Method:</b> Protein A purification
<b>Source:</b> Rabbit	<b>GeneID (NCBI):</b> 4082	<b>CloneNo.:</b> 253128F1
<b>Isotype:</b> IgG	<b>UNIPROT ID:</b> P29966	
<b>Immunogen Catalog Number:</b> EG6453	<b>Full Name:</b> myristoylated alanine-rich protein kinase C substrate	
	<b>Calculated MW:</b> 32 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



The mean MARCKS concentration was determined to be 83.8 ng/mL in HEK-293 cell extract based on a 1.2 mg/mL extract load, 32.3 ng/mL in SH-SY5Y cell extract based on a 1.5 mg/mL extract load, 8.3 ng/mL in A-204 cell extract based on a 1.2 mg/mL extract load.



Cytometric bead array standard curve of MP03172-1, MARCKS Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87693-3-PBS. Detection antibody: 87693-2-PBS. Standard: Eg6453. Range: 0.312-40 ng/mL.