

For Research Use Only

MARCS Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 87693-3-PBS



Basic Information

Catalog Number:

87693-3-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG6453

GenBank Accession Number:

NM_002356.7

GeneID (NCBI):

4082

UNIPROT ID:

P29966

Full Name:

myristoylated alanine-rich protein kinase C substrate

Calculated MW:

32 kDa

Purification Method:

Protein A purification

CloneNo.:

253128D11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean MARCKS concentration was determined to be 83.8 ng/mL in HEK-293 cell extract based on a 1.2 mg/mL extract load, 32.3 ng/mL in SH-SY5Y cell extract based on a 1.5 mg/mL extract load, 8.3 ng/mL in A-204 cell extract based on a 1.2 mg/mL extract load.



Cytometric bead array standard curve of MP03172-1, MARCKS Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87693-3-PBS. Detection antibody: 87693-2-PBS. Standard: Eg6453. Range: 0.312-40 ng/mL.