

For Research Use Only

ABCC11 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 87699-1-PBS



Basic Information

Catalog Number: 87699-1-PBS	GenBank Accession Number: BC039085	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 85320	CloneNo.: 253118C9
Isotype: IgG	UNIPROT ID: Q96J66	
	Full Name: ATP-binding cassette, sub-family C (CFTR/MRP), member 11	
	Calculated MW: 150 kDa, 120 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean ABCC11/MRP8 concentration was determined to be 3.5 ng/mL in MCF-7 cell extract based on a 1.2 mg/mL extract load, U-87 MG cell extract were served as a negative control.



Cytometric bead array standard curve of MP03159-1, ABCC11/MRP8 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87699-2-PBS. Detection antibody: 87699-1-PBS. Standard: SY03027. Range: 0.078-20 ng/mL