

For Research Use Only

ABCC11 Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 87699-2-PBS



Basic Information

Catalog Number:

87699-2-PBS

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC039085

GeneID (NCBI):

85320

UNIPROT ID:

Q96J66

Full Name:

ATP-binding cassette, sub-family C (CFTR/MRP), member 11

Calculated MW:

150 kDa, 120 kDa

Purification Method:

Protein A purification

CloneNo.:

253118G8

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean ABCC11/MRP8 concentration was determined to be 3.5 ng/mL in MCF-7 cell extract based on a 1.2 mg/mL extract load, U-87 MG cell extract were served as a negative control.



Cytometric bead array standard curve of MP03159-1, ABCC11/MRP8 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87699-2-PBS. Detection antibody: 87699-1-PBS. Standard: SY03027. Range: 0.078-20 ng/mL