

For Research Use Only

# IL17C Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number:87706-1-PBS



## Basic Information

Catalog Number:

87706-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG2788

GenBank Accession Number:

NM\_013278.3

GeneID (NCBI):

27189

UNIPROT ID:

Q9P0M4

Full Name:

interleukin 17C

Calculated MW:

22 kDa

Purification Method:

Protein A purification

CloneNo.:

251249F9

## Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



HCT-116 cells were cultured unstimulated or stimulated with 15 ng/mL IL-1 $\beta$  for 1 day. The mean IL-17C concentration was determined to be 10.2 pg/mL in unstimulated cell culture supernatant, 453.2 pg/mL in stimulated cell culture supernatant.



Sandwich ELISA standard curve of MP03144-1, Human IL-17C Recombinant Matched Antibody Pair - PBS only. 87706-2-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg2788. 87706-1-PBS was HRP conjugated as the detection antibody. Range: 15.6-1000 pg/mL