

For Research Use Only

Mouse TREML4/TLT-4 Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 87709-2-PBS



Basic Information

Catalog Number: 87709-2-PBS	GenBank Accession Number: NM_172623.2	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 224840	CloneNo.: 252981D5
Isotype: IgG	UNIPROT ID: Q3LRV9-1	
Immunogen Catalog Number: EG7514	Full Name: triggering receptor expressed on myeloid cells-like 4	
	Calculated MW: 29 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
mouse

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean TLT-4/TREML4 concentration was determined to be 98.2 ng/mL in mouse spleen tissue extract based on a 6.4 mg/mL extract load, 2.9 ng/mL in mouse lung tissue extract based on a 5.0 mg/mL extract load.



Cytometric bead array standard curve of MP03163-1, MOUSE TLT-4/TREML4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87709-2-PBS. Detection antibody: 87709-1-PBS. Standard: Eg7514. Range: 0.781-25 ng/mL.