

For Research Use Only

# FOLR3 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number:87719-1-PBS



## Basic Information

Catalog Number: 87719-1-PBS	GenBank Accession Number: NM_000804.3	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 2352	CloneNo.: 250411D7
Isotype: IgG	UNIPROT ID: P41439	
Immunogen Catalog Number: EG3170	Full Name: folate receptor 3 (gamma)	
	Calculated MW: 28kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



The mean FOLR3 concentration was determined to be 1.5 ng/mL in human serum 1, 3.5 ng/mL in human serum 2, 15.7 ng/mL in human serum 3.



Cytometric bead array standard curve of MP03170-1, FOLR3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87719-2-PBS. Detection antibody: 87719-1-PBS. Standard: Eg3170. Range: 0.156-20 ng/mL.