

For Research Use Only

# HRH1 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 87723-1-PBS



## Basic Information

Catalog Number:

87723-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG30464

GenBank Accession Number:

BC060802

GeneID (NCBI):

3269

UNIPROT ID:

P35367

Full Name:

histamine receptor H1

Calculated MW:

487 aa, 56 kDa

Purification Method:

Protein A purification

CloneNo.:

253277E3

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



The mean HRH1 concentration was determined to be 12.9 ng/mL in U-87MG cell extract based on a 1.2 mg/mL extract load, 19.5 ng/mL in Jurkat cell extract based on a 1.2 mg/mL extract load, 3.0 ng/mL in Ramos cell extract based on a 1.5 mg/mL extract load.



Cytometric bead array standard curve of MP03173-1, HRH1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87723-2-PBS. Detection antibody: 87723-1-PBS. Standard: Ag30464. Range: 0.625-80 ng/mL.