

For Research Use Only

CCDC80 Recombinant monoclonal antibody, PBS Only

Catalog Number: 87731-2-PBS



Basic Information

Catalog Number: 87731-2-PBS	GenBank Accession Number: BC116178	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 151887	CloneNo.: 253134D6
Isotype: IgG	UNIPROT ID: Q76M96	
Immunogen Catalog Number: EG6754	Full Name: coiled-coil domain containing 80	
	Calculated MW: 950 aa, 108 kDa	
	Observed MW: 110-180 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human

Background Information

Coiled-coil domain-containing protein 80 (CCDC80) is a member of the coiled-coil domain-containing protein family, which is mainly expressed in vascular smooth muscle cells. It plays an integral role in the regulation of cardiovascular remodeling and homeostasis (PMID: 23831687; 30439364; 27375481). In SDS PAGE, CCDC80 can be detected as double bands (PMID: 40278823; 35843176).

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

