

For Research Use Only

Mouse CD222/Igf2r Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 87733-1-PBS



Basic Information

Catalog Number:

87733-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG6300

GenBank Accession Number:

NM_010515.2

GeneID (NCBI):

16004

UNIPROT ID:

Q07113

Full Name:

insulin-like growth factor 2 receptor

Calculated MW:

274 kDa

Purification Method:

Protein A purification

CloneNo.:

253081C9

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean IGF-II R/IGF2R concentration was determined to be 565.1 ng/mL in NIH/3T3 cell extract based on a 1.5 mg/mL extract load, 769.7 ng/mL in mouse lung tissue extract based on a 5.0 mg/mL extract load.



Cytometric bead array standard curve of MP03179-1, MOUSE IGF-II R/IGF2R Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87733-2-PBS. Detection antibody: 87733-1-PBS. Standard: Eg6300. Range: 0.625-80 ng/mL.