

For Research Use Only

MLKL Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 87741-2-PBS



Basic Information

Catalog Number: 87741-2-PBS	GenBank Accession Number: BC028141	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 197259	CloneNo.: 253011H6
Isotype: IgG	UNIPROT ID: Q8NB16	
Immunogen Catalog Number: AG15188	Full Name: mixed lineage kinase domain-like	
	Calculated MW: 471 aa, 54 kDa	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



The mean MLKL concentration was determined to be 6.73 ng/mL in HUVEC cell extract based on a 1.2 mg/mL extract load.



Plasma of six individual healthy human donors was measured. The human MLKL concentration of detected samples was determined to be 2.02 ng/mL with a range of 0.22 - 6.91 ng/mL.



Sandwich ELISA standard curve of MP03218-1, Human MLKL Recombinant Matched Antibody Pair - PBS only. 87741-1-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag15188. 87741-2-PBS was HRP conjugated as the detection antibody. Range: 0.195-12.5 ng/mL.