

For Research Use Only

# PYK2 Recombinant monoclonal antibody, PBS Only

Catalog Number: 87746-1-PBS



## Basic Information

Catalog Number:

87746-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG11487

GenBank Accession Number:

BC042599

GeneID (NCBI):

2185

UNIPROT ID:

Q14289

Full Name:

PTK2B protein tyrosine kinase 2 beta

Calculated MW:

1009 aa, 116 kDa

Observed MW:

112-115 kDa

Purification Method:

Protein A purification

CloneNo.:

253299F11

## Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse, rat

## Background Information

Proline-rich tyrosine kinase 2 (Pyk2; also known as CAK, RAFTK and CADTK) is a cytoplasmic tyrosine kinase implicated to play a role in several intracellular signaling pathways. It is expressed in many cells and tissues and migrated as a 130 kDa band in Western Blotting analysis. The size of PYK2 is 115 kDa, and it expressed in many cells and tissues and migrated as a 130 kDa band, but several lower molecular weight bands (75 kDa, 80 kDa, and 97 kDa) were seen in Western Blotting analysis, possibly due to proteolytic degradation.

## Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 87746-1-RR (PYK2 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 87746-1-PBS in a different storage buffer formulation.