

For Research Use Only

ANXA6 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 87752-2-PBS



Basic Information

Catalog Number:

87752-2-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3263

GenBank Accession Number:

BC017046

GeneID (NCBI):

309

UNIPROT ID:

P08133

Full Name:

annexin A6

Calculated MW:

673 aa, 76 kDa

Purification Method:

Protein A purification

CloneNo.:

253305E9

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean Annexin VI concentration was determined to be 675.5 ng/mL in HeLa cell extract based on a 1.2 mg/mL extract load, 1316.3 ng/mL in A-673 cell extract based on a 1.5 mg/mL extract load, HaCaT cell extract were served as a negative control.



Cytometric bead array standard curve of MP03192-1, Annexin VI Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87752-3-PBS. Detection antibody: 87752-2-PBS. Standard: Ag3263. Range: 1.562-200 ng/mL.