

For Research Use Only

Cystatin F Recombinant monoclonal antibody, PBS Only

Catalog Number:87771-1-PBS



Basic Information

Catalog Number: 87771-1-PBS	GenBank Accession Number: BC015507	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 8530	CloneNo.: 253239D1
Isotype: IgG	UNIPROT ID: O76096	
Immunogen Catalog Number: EG8049	Full Name: cystatin F (leukocystatin)	
	Calculated MW: 167 aa, 19 kDa	
	Observed MW: 16-25 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human

Background Information

Cystatin F (also known as CST7, CMAP or leucocyte cystatin) is an endogenous cysteine protease inhibitor belonging to type II cystatin superfamily. Covalently linked homodimer is formed through C26-C63 disulfide bonds, which is the main form of its secretion. It contains N- glycosylation sites (Asn62 and Asn115), and glycosylation modification determines its intracellular localization-high mannose type guides it to target lysosomes, and complex sugar chain type promotes secretion. The active form is monomer, and the dimer needs to be activated in lysosomal reduction environment, and the 15 amino acids at the N-terminal are removed by cathepsin V to become a completely active form. Mainly expressed in peripheral blood cells (such as cytotoxic lymphocytes, eosinophils, neutrophils, mast cells) and spleen.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 87771-1-RR (Cystatin F antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 87771-1-PBS in a different storage buffer formulation.