

For Research Use Only

Mouse ORM1 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 87784-1-PBS



Basic Information

Catalog Number:

87784-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG2752

GenBank Accession Number:

NM_008768.2

GeneID (NCBI):

18405

UNIPROT ID:

Q60590

Full Name:

orosomuroid 1

Calculated MW:

24 kDa

Purification Method:

Protein A purification

CloneNo.:

253420F8

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean ORM1 concentration was determined to be 27.4 ug/mL in mouse serum 1, 21.4 ug/mL in mouse serum 2.



Cytometric bead array standard curve of MP03203-1, MOUSE ORM1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87784-2-PBS. Detection antibody: 87784-1-PBS. Standard: Eg2752. Range: 0.625-80 ng/mL.