

For Research Use Only

EDEM1 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 87789-2-PBS



Basic Information

Catalog Number: 87789-2-PBS	GenBank Accession Number: BC019088	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 9695	CloneNo.: 253376A6
Isotype: IgG	UNIPROT ID: Q92611	
Immunogen Catalog Number: AG41790	Full Name: ER degradation enhancer, mannosidase alpha-like 1	
	Calculated MW: 74 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean EDEM1 concentration was determined to be 18.0 ng/mL in L02 cell extract based on a 1.2 mg/mL extract load. Caco-2 cell extract were served as a negative control.



Cytometric bead array standard curve of MP03204-1, EDEM1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87789-1-PBS. Detection antibody: 87789-2-PBS. Standard: Ag41790. Range: 0.469-60 ng/mL.