

For Research Use Only

MARCH5 Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 87790-3-PBS



Basic Information

Catalog Number: 87790-3-PBS	GenBank Accession Number: BC015480	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 54708	CloneNo.: 253335F9
Isotype: IgG	UNIPROT ID: Q9NX47	
Immunogen Catalog Number: AG2806	Full Name: membrane-associated ring finger (C3HC4) 5	
	Calculated MW: 278 aa, 31 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean MARCH5 concentration was determined to be 9.6 ng/mL in A431 cell extract based on a 1.5 mg/mL extract load, 4.9 ng/mL in THP-1 cell extract based on a 1.2 mg/mL extract load.



Cytometric bead array standard curve of MP03206-1, MARCH5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87790-3-PBS. Detection antibody: 87790-2-PBS. Standard: Ag2806. Range: 0.625-80 ng/mL.