

For Research Use Only

C1QL2 Recombinant monoclonal antibody, PBS Only

Catalog Number:87797-1-PBS



Basic Information

Catalog Number: 87797-1-PBS	GenBank Accession Number: NM_182528	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 165257	CloneNo.: 253329D5
Isotype: IgG	UNIPROT ID: Q7Z5L3	
Immunogen Catalog Number: AG37614	Full Name: complement component 1, q subcomponent-like 2	
	Calculated MW: 29 kDa	
	Observed MW: 35 kDa	

Applications

Tested Applications:
WB, IF/ICC, Indirect ELISA

Species Specificity:
human, mouse, rat, pig

Background Information

Complement C1q-like protein 2 (C1QL2, also known as CTRP10) is a member of the C1q-like protein family. It is a secreted protein that plays a role in synaptic organization and function. C1QL2 is involved in the formation and maintenance of synapses, particularly in the central nervous system (CNS) (PMID: 20646056). It interacts with specific neurexin-3 isoforms to recruit synaptic vesicles and modulate synaptic transmission. Additionally, C1QL2 has been implicated in the regulation of hippocampal mossy fiber-CA3 synapse function, suggesting its role in synaptic plasticity and memory formation (PMID: 28219683).

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 87797-1-RR (C1QL2 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 87797-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (-20°C Methanol) fixed HeLa cells using C1QL2 antibody (87797-1-RR, Clone: 253329D5) at dilution of 1:400 and CoraLite® 488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2). This data was developed using the same antibody clone with 87797-1-PBS in a different storage buffer formulation.