

For Research Use Only

ISG15 Recombinant monoclonal antibody, PBS Only

Catalog Number: 87799-2-PBS



Basic Information

Catalog Number: 87799-2-PBS	GenBank Accession Number: BC009507	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 9636	CloneNo.: 260123E6
Isotype: IgG	UNIPROT ID: P05161	
Immunogen Catalog Number: AG8782	Full Name: ISG15 ubiquitin-like modifier	
	Calculated MW: 165 aa, 18 kDa	
	Observed MW: 15 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human

Background Information

ISG15 is a ubiquitin-like protein that becomes conjugated to many cellular proteins upon activation by interferon- α (IFN α) and - β (IFN β). ISG15 forms covalent conjugates with its target proteins in a process called ISGylation, which in mammals is known to play a role in antiviral immunity. ISG15 proteins possess two ubiquitin-like (UBL) domains and a highly conserved C-terminal LRG sequence, the latter being known as the ubiquitin conjugation motif. Intracellular ISG15 are conjugated, via the LRG motif, to target proteins through a process called ISGylation, which resembles largely ubiquitination, the process of formation of ubiquitin conjugates. Unconjugated extracellular ISG15, which are released from several types of human and murine cells, are known to possess cytokine-like activity.

Storage


Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 87799-2-RR (ISG15 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 87799-2-PBS in a different storage buffer formulation.