

For Research Use Only

FRS3 Recombinant monoclonal antibody

Catalog Number: 87810-1-RR



Basic Information

Catalog Number:

87810-1-RR

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2896

GenBank Accession Number:

BC010611

GeneID (NCBI):

10817

UNIPROT ID:

O43559

Full Name:

fibroblast growth factor receptor substrate 3

Calculated MW:

492 aa, 54 kDa

Observed MW:

54 kDa

Purification Method:

Protein A purification

CloneNo.:

253396B5

Recommended Dilutions:

WB: 1:5000-1:50000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : DU 145 cells, mouse brain tissue, PC-3 cells, rat brain tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 87810-1-RR (FRS3 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.



Various lysates were subjected to SDS PAGE followed by western blot with 87810-1-RR (FRS3 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.